

Style Follow Skinable Face Cooler

Cooling Massager That Soothes Puffiness and Contracts Pores on Your Skin



- Boosting Elasticity



- Cooling Therapy
- Pore Tightening Reduce Skin Temperature



- Absorption of Cosmetics

Skinable Face cooler Tips



- Accessory Holder with High Utilization
- Chilliness Lasts about 1 Hour (* Depending on the room temp.)
- Standing Out Cute Design
- Compact Size to Carry Anywhere

· HOW TO USE ·

- 1. Separate refrigerant part from the main body, then keep it in the freezer over 1 hour.
- 2. Assemble it with the main body then gently massage on your face or body.
 - * You can use with soothing gel or cream for better effectiveness.

SBS Style Follow Official Sponsorship

Han TheHαn Co., Ltd. Website http://sf.skinable.net Tel. 032.579.3833 Email thehan@thehan.net