

Real Multi Beauty Device Skingble True Solution LED Massager

Galvanic Ion Therapy

What's Galvanic?

Galvanic technology works by using positive and negative currents. The benefits include smoother and clearer skin, as well as reduced wrinkles and skin tightening.

Deep Cleansing (Positive(+) Ion)

Positive electrons attract dirt and impurities from the deep layers and pores to the skin's surface.

Nutrition (Negative(-) Ion)

The negative electrons ensure the absorption of nutrients deep into the skin's lavers.



Red (625nm) : Promotion of collagen production,

Boosting elasticity

Blue (465nm): Skin health improvements.

Remove bactreia causina acne

Green (525nm): Skin tone improvements. Calming irritated skin

Microwave Soothing Massage

Micro-wave massage mode vibrates 12,000 times per minute. It promotes skin elasticity and V-line by massaging tight muscles of face,

Ultrasound Massage (3Mhz)

Ultrasound delivers heat energy to deep in the dermal layer. causing a small injury. The body responds to this injury by creating more collagen, which in turn helps lift and tighten the skin.





LED Therapy



Tightening



Brightening

